BioGenex

8810 Kato Road, Suite 100E & 200E,
Fremont,CA94538

Tel : +1(800)421-4149,
Fax: +1(510)824-1490,
support@biogenex.com

EC REP

EmergoEurope,
Prinsessegracht
20 2514 AP

The Hague,
The Netherlands

BLOCKING REAGENTS

Doc. No. 932-HK083E-4 Rev. No. F
Release Date: 10-Aug-2020

ENGLISH

Intended Use

BioGenex Blocking reagents are intended for in vitro diagnostic use. They are used for inhibition and/or quenching of non-specific targets in formalin-fixed tissue sections in immunohistochemistry (IHC) or in situ hybridization (ISH) assays.

Summary and Explanation

Elimination or reduction of non-specific background staining is an essential step in IHC or ISH procedures. BioGenex offers a variety of blocking agents such as power block, protein block, peroxide block, avidin-biotin block, blocking agents I and II and levamisole, each intended for inhibiting a specific group of non-specific targets.

Principles of the Procedure

Blocking reagents bind to or saturate non-specific targets in formalin-fixed paraffin embedded tissues thereby preventing non-specific background staining. This process helps in enhancing specific antigen antibody reactions in IHC or ISH procedures. The blocking reagents and their non-specific targets are shown in the following table:

Blocking reagent	Non-specific targets
Power Block	Saturates non specific positive charges both on slide as well as in tissue sections
Peroxide Block	Saturates the endogenous peroxidases in tissue sections
Avidin-Biotin blocking kit	Binds to endogenous biotin, biotin receptors or avidin binding sites in tissue sections
Protein block / normal goat serum	Saturates the non-specific goat immunoglobulin binding sites in tissue sections
Protein block / normal rabbit serum	Saturates the non-specific rabbit immunoglobulin binding sites in tissue sections
Levamisole	Competitive inhibition of endogenous alkaline phosphatase in tissue sections

Reagents Provided

Catalog No.	Blocking reagent
HK085-5KE	Power Block
HK111-5K HK111-XAK	Peroxide Block
HK102-10KE HK102-5KE HK102-20XE	Avidin-Biotin blocking kit
HK112-9KE	Protein block / normal goat serum

Reagents Required but Not Supplied

All reagents required for blocking step are provided. See antibody and detection kit datasheet for complete set of reagents required for immunohistochemistry (IHC) or in situ hybridization (ISH) procedures.

Storage and Handling

Store all reagents at 2-8°C. Do not use after expiration dates indicated on the reagent labels.

Specimen Preparation

Paraffin embedded and formalin fixed tissue sections 4-6 micron thick are suitable for use with enzyme pretreatment.

Precautions

Power Block™, Protein Block and Avidin-Biotin Blocking Kit contain sodium azide, at concentrations of less than 0.1%. Sodium azide is not classified as a hazardous chemical at the product concentrations. However, toxicity information regarding sodium azide at product concentrations has not been thoroughly investigated. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing (Center for Disease Control, 1976, National Institute

for Occupational Safety and Health, 1976). For more information, a Material Safety Data Sheet for sodium azide in pure form is available upon request.

Power Block™ is not appropriate as a wash buffer in procedure steps immediately prior to application of horseradish peroxidase (HRP) label.

Power Block™ is a potent blocker and it should not remain on the tissue for more than ten minutes.

Preparation of working solution

Levamisole: Prepare a fresh solution by adding 100 micro liter of 50X concentrated solution to 5 ml of alkaline phosphatase substrate solutions.

Power Block, 10X: Dilute 10 fold using distilled water.

Staining Protocol

The following table provides the recommended treatments with blocking reagents. Refer to the antibody and detection kit datasheets for complete staining protocols. User should validate any other conditions. For automated follow instructions in operator's manual.

Blocking reagent	BioGenex Recommendations
Power Block	Incubate for5-10 min at room temperature. No wash step required prior to incubation with antibody.
Peroxide Block	Incubate for 10 min at room temperature after dewaxing and antigen retrieval and wash in appropriate buffer
Avidin-Biotin blocking kit	Incubate avidin block for 15 min, wash and follow by biotin block for 15 min at room temperature and wash.
Protein block / normal goat serum	Incubate for 20 min at room temperature followed by incubation with primary antibody. No washing required prior to incubation with antibody.
Protein block / normal rabbit serum	Incubate for 20 min at room temperature followed by incubation with primary antibody. No washing required prior to incubation with antibody.
Levamisole	Prepare a fresh solution by adding 100 micro liter 50X concentrated solution to 5ml of alkaline phosphatase. substrate solutions. Incubate on slide following instructions for the substrate solution.

Quality Control

Refer to the appropriate detection system package inserts for guidance on general quality control procedures.

Trouble shooting

Run appropriate positive and negative controls to ensure the optimal signal.

Refer to the troubleshooting section in the package inserts of BioGenex Super Sensitive Detection Systems (or other equivalent detection systems) for remedial actions on detection system related issues, or contact BioGenex Technical Service Department at (925) 275-0550 to report unusual staining.

Expected Results

The non-specific background in the IHC or ISH stained slides would be is either completely eliminated or negligible.

The colored pigment deposition will be seen at the specific target site. The color signal can be visualized using light microscopy.

Limitations of the procedure

Incomplete blocking may lead to background staining. Incubation of power block longer than 10 minutes may lead to reduced signal intensity.

Performance Characteristics

BioGenex has conducted studies to evaluate the performance of all its blocking reagents using several BioGenex antibodies and detection systems. The blocking agents have been found show minimal to no non-specific binding thereby enhancing specific antigen antibody reactions. BioGenex blocking reagents have shown reproducible and consistent results when used within a single run, between runs, between lots and wherever applicable between manual and automated runs. The products have been determined to be stable for the periods specified on the labels either by standard real time or accelerated testing methods. BioGenex ensures product quality through 100% quality control for all products released and through surveillance programs.

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Bibliography / Bibliografía / Bibliografie / Bibliografia

- SheehanDC and Hrapchak BB. Theory and Practice of Histotechnology. St. Louis: C.V. Mosby Co. 1980.
- Kiernan JA. Histological and Histochemical Methods: Theory and Practice. New York: Pergamon Press 1981.
- Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. Lab Med 1983; 14:767-770.
- Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immuno histochemistry. Am J ClinPathol 1980 May;73(5):626-632.
- National Institute for Occupational Safety and Health,(NIOSH), Rockville, MD. Explosive azide hazard, Publication No. 78-127, 1976.
- Center for Disease Control. Decontamination of Laboratory Sink Drains to Remove Azide Salts. Center for Disease Control Manual Guide--Safety Management, No. CDC-22, Atlanta, Georgia. April 30, 1976.

EC|REP Representative in the European Community
Mandatario nella Comunità Europea
Bevollmächtigter in der Europäischen Gemeinschaft
Representante autorizado en la Comunidad Europea

IVD

In Vitro Diagnostic Medical Device
Dispositivo medico-diagnostico in vitro
In Vitro Diagnostikum
Producto sanitario para diagnóstico in vitro

i

Cons ult Instructions for use
Consultare le istruzioni per l'uso
Gebrauchsanweisungbeachten
Consulte las instruccionesde uso

T

Temperature Limitation
Limiti di temperatura
ZulässigerTemperaturbereich
Limite de temperatura

M

Manufacturer
Fabbriante
Hersteller
Fabricante