








## Super Sensitive™ Polymer-HRP IHC Detection For Manual and Automated Systems

### Kit Components:

Cat. No.	No. of Slides	EZ-AR™ 1	EZ-AR™ 2	Peroxid eBlock	Power Block	Mouse Negative Control	Rabbit Negative Control	Super Enhancer™	Poly-HRP	Substrate	Chromogen	Empty Mixing Vial	Hematoxylin	Xmount	Symbol
<b>QD400-60KE</b> Manual	<b>60</b>	---	---	HK111-5KT (1x6 ml)	HK083-5KT (1x6 ml)	HK119-5M (1x3 ml)	HK408-5RT (1x3ml)	HK518-06K (1x6 ml)	HK519-06K (1x6 ml)	HK520-10K (1x10 ml)	HK124-1K (1x1 ml)	Mixing Vial (1x6ml)	HK100-5KT (1x6 ml)	---	
<b>QD420-YIKE</b> Manual	<b>500</b>	---	---	---	---	---	---	HK518-50K (1x50 ml)	HK519-50K (1x50 ml)	HK520-50K (1x50 ml)	HK124-025K (1 x 2.5 ml)	Mixing Vial (1x6ml)	---	---	
<b>QD430-XAKE</b> Manual	<b>1000</b>	---	---	---	---	---	---	HK518-YAK (1x100 ml)	HK519-YAK (1x100 ml)	HK520-YAK (1x100 ml)	HK124-05K (1x5 ml)	Mixing Vial (1x6ml)	---	---	
<b>QD440-XAKE</b> Manual	<b>1000</b>	---	---	---	---	---	---	HK518-YAK (1x100 ml)	HK519-YAK (1x100 ml)	---	---	---	---	---	
<b>QD550-YCXE</b> Xmatrix Infinity	<b>200</b>	HX031-YCX (1x16 ml)	HX032-YCX (1x16ml)	HX026-YAX (3x16 ml)	HX083-YAX (3x 16 ml)	---	---	HX027-YCX (2x16 ml)	HX028-YCX (2x16 ml)	HX029-50X (4 x 13 ml)	HX010-YCK (1x4ml)	HX616-YCX (4no)	HX030-YAX (3 x 16 ml)	HX035-10XT (1x15ml)	
<b>QD550-YCDE</b> Xmatrix Elite	<b>200</b>	HX031-YCDE (1x16 ml)	HX032-YCDE (1x16 ml)	HX026-YADE (3 x 16 ml)	HX083-YADE (3 x 16 ml)	---	---	HX027-YCDE (2 x 16 ml)	HX028-YCDE (2 x 16 ml)	HX029-50DE (4 x 13 ml)	HX010-07XE (1x4ml)	HX615-YADE (4no)	HX030-YADE (3 x 16 ml)	HX035-YCDE (1x15ml)	
<b>QD410-YAXE</b> i6000 Diagnostics	<b>200</b>	---	---	HK111-20K (20 ml)	HK083-20K (20 ml)	---	---	HK518-20X (20 ml)	HK519-20X (20 ml)	HK520-10X (5 x 10ml)	HK124-7KT (1 x 4ml)	---	HK100-20K (20ml)	---	

**Note: It is recommended that reagents not be substituted across kit lot numbers.**

©2020, BioGenex Laboratories. All rights reserved.

<b>Category</b>	Detection Systems	<b>Revision No.</b>	C
<b>Document No.</b>	932-IHC2ST-Appendix	<b>Release Date</b>	04-Jan-2022