






SUPER SENSITIVE™* ISH DETECTION SYSTEM (One Step)

Product Name	Cat. No.	No. of Slides	Liquid Pepsin	NAR	Hybridization Solution II	Wash Buffer A	Wash Buffer B	Wash Buffer E	Wash Buffer F	Peroxi de Block	Power Block	Anti-Fluorescent Antibody	One Step Poly-HRP	Substrate	Chromogen	Hemat oxylin	Mixing Vial Empty	Symbol
ISH Detection Kit- DEMO	DF400-25K	25	HK632-01K (2.5 ml)	HK601-01K (1.5ml)	HK881-3K (3 ml)	HK839-06K (10 ml)	HK880-06K (10 ml)	HK946-06K (10 ml)	HK947-06K (10 ml)	HK026-03K (3 ml)	HK083-03K (3 ml)	HK818-01K (1.5 ml)	HK943-01K (1.5 ml)	HK520-06K (5 ml)	HK124-005K (0.5 ml)	HK030-03K (3 ml)	Mixing Vial (6 ml)	
XISH™ One Step Polymer HRP Detection System for Xmatrx	DF400-50KE	50	HK632-05KE (5 ml)	HK601-03KE (3 ml)	HK881-5KE (5 ml)	HK839-20KE (20 ml)	HK880-20KE (20 ml)	HK946-20KE (20 ml)	HK947-20KE (20 ml)	HK026-05KE (5.5ml)	HK083-5KE (5.5ml)	HK818-03KE (3 ml)	HK943-03KE (3 ml)	HK520-10KE (10 ml)	HK124-15KT (1.5 ml)	HK030-10KE (5 ml)	Mixing Vial (6 ml)	
XISH™ One Step Polymer HRP Detection System for Xmatrx	DF400-YADE	100	HK632-04XE (5 ml)	HX601-5XE (5 ml)	HK881-5KE (6 ml)	HX839-10XE (10 ml x 2)	HX880-10XE (10 ml x 2)	HX946-10XE (10 ml x 2)	HX947-10XE (10 ml x 2)	HX026-10XE (10 ml)	HX083-10XE (10 ml)	HX818-04XE (5 ml)	HX943-04XE (5 ml)	HX029-05XE (5 ml X 4)	HX010-02XE (2ml)	HK030-10XE (10 ml)	HX615-YADE (16 ml X 5)	

*SS=Super Sensitive

SUPER SENSITIVE™* ISH DETECTION SYSTEM (One Step)

XISH™ One step Polymer- HRP Detectio n Kit for NanoVIP (Open System)	DF540 -YADX	100	HX051- 06X	HX601- 08X	HX881- 08X	HX839- 06X	HX880- 06X	HX946- 06X	HX947- 06X	HX026- 06X	HX083- 06X	HX818-08x	HX943- 06X	HX029- 08X	HX124- 04X	HX030- 06X	HX616- 10X	
			(6ml x2)	(8ml x1)	(8ml x1)	(6ml x2)	(6ml x2)	(6ml x2)	(6ml x2)	(6ml x2)	(6ml x2)	(6ml x2)	(8ml x1)	(6ml x2)	(8ml x2)	(4ml x1)	(6ml x2)	
XISH™ One step Polymer- HRP Detectio n Kit for NanoVIP (Close System)	DF541 -YADX	100	HX041- 06X	HX606- 08X	HX885- 08X	HX838- 06X	HX886- 06X	HX941- 06X	HX944- 06X	HX037- 06X	HX042- 06X	HX819-08X	HX942- 06X	HX038- 08X	HX125- 04X	HX039- 06X	HX617- 10X	
			(6ml x2)	(8ml x1)	(8ml x1)	(6ml x2)	(6ml x2)	(6ml x2x)	(6ml x2)	(6ml x2)	(6ml x2)	(6ml x2)	(8ml x1)	(6ml x2)	(8ml x2)	(4ml x1)	(6ml x2)	

*SS=Super Sensitive