

DATA SHEET

Xmatrx™ Coverslip

Doc. No.932-XT118, Rev. No. C
Date of Release: 15-Sep-2021

Intended Use

This product is intended for laboratory use only. BioGenex Coverslips are to be used on Xmatrx® Staining system.

Specifications

1. Coverslips are used on Xmatrx® Staining system after every time reagent is applied to optimize use and minimize wastage of micro-reagents.
2. Coverslips come in three configurations to accommodate to fit with different Micro Chamber slides provided by Biogenex: a single full-size test area of 25 x 40 mm, a single 2/3-size test area of 25 x 25 mm, and three 1/3-size test areas per slide, each compartment measuring 18 x 18 mm.

Items Supplied

One of the following:

| Cat. Nos. | Description |
|------------------|------------------------------|
| XT118-50X | 25 x 40 Coverslips (50 Nos) |
| XT118-YRK | 25 x 40 Coverslips (500 Nos) |

Applications

1. Coverslips are to be used with Xmatrx® Automated Staining system.
2. Coverslip application optimizes use and minimizes wastage of micro-reagents.
3. It also avoids the evaporation of micro-reagents during incubation period.