

## DATA SHEET

### *Xmatrx® Coverslip*

---

Doc. No.932-XT121, Rev. No. A  
Date of Release: 13-May-2014

#### **Intended Use**

This product is intended for laboratory use only. BioGenex Coverslips are to be used on Xmatrx® Staining system.

#### **Specifications**

1. Coverslips are used on Xmatrx® Staining system after every time reagent is applied to optimize use and minimize wastage of micro-reagents.
2. Coverslips come in three configurations to accommodate to fit with different barrier slides provided by Biogenex: a single full-size test area of 25 x 40 mm, a single 2/3-size test area of 25 x 25 mm, and three 1/3-size test areas per slide, each compartment measuring 18 x 18 mm.

#### **Items Supplied**

One of the following:

<b>Cat. Nos.</b>	<b>Description</b>
XT121-YBX	18x18 mm, Coverslips for Xmatrx® (175 Coverslips/Box)
XT121-XBX	18x18 mm, Coverslips for Xmatrx® (1750 Coverslips/Case)

#### **Applications**

1. Coverslips are to be used with Xmatrx® Automated Staining system.
2. Coverslip application optimizes use and minimizes wastage of micro-reagents.
3. It also avoids the evaporation of micro-reagents during incubation period.