

DATA SHEET

Xmatrx® Infinity Coverslip

Doc. No.932-XT122, Rev. No. A
Date of Release: 13-May-2014

Intended Use

This product is intended for laboratory use only. BioGenex Coverslips are to be used on Xmatrx® Infinity Staining system.

Specifications

1. Coverslips are used on Xmatrx® Staining Infinity system after every time reagent is applied to optimize use and minimize wastage of micro-reagents.
2. Coverslips come in three configurations to accommodate to fit with different barrier slides provided by Biogenex: a single full-size test area of 25 x 40 mm, a single 2/3-size test area of 25 x 25 mm, and three 1/3-size test areas per slide, each compartment measuring 18 x 18 mm.

Items Supplied

One of the following:

Cat. Nos.	Description
XT122-90X	Coverslips, 25x25 mm, Xmatrx® Infinity 90 Coverslips/Box)
XT122-YQK	Coverslips, 25x25 mm, Xmatrx® Infinity(500 Nos)

Applications

1. Coverslips are to be used with Xmatrx® Automated Staining system.
2. Coverslip application optimizes use and minimizes wastage of micro-reagents.
3. It also avoids the evaporation of micro-reagents during incubation period.