

DATA SHEET**Xmatrix, i6000™ and OptiMax® Pipette Tips**

Doc. No.932-XT145, Rev. No: D

Date of Release: 25-Oct-2019

Intended Use

This product is intended for laboratory use only. BioGenex pipette tips are optimized for use on BioGenex Xmatrix®, i6000™ and OptiMax® Staining Systems.

Specifications

BioGenex pipette tips are made of high-quality polypropylene and are RNase and heavy metals-free when untampered. Inner surface is extremely smooth and requires minimum wetting.

Items Supplied

One of the following:

Cat. No.	Description
XT145-05X	1 case (4800 tips) of 200 µl Xmatrix® FISH pipette tip
XT146-01X	1 box (960 tips) of 200 µl Xmatrix® FISH pipette tip
XT116-XAX	1 case (960 tips) of 200 µl Pipette tips Xmatrix®
XT116-PT	1 case (96 tips) of 200 µl Pipette tips Xmatrix®
XT104-05X	1 case (960 tips) of 1ml Pipette tips i6000™
XT105-01X	1 box (192 tips) of 1ml Pipette tips for i6000™
XT106-PT	1 case (1000 tips) of 1ml Pipette tips for OptiMax® Upgrade
XT107-PT	1 box (100 tips) of 1ml Pipette tips for OptiMax® Upgrade
XT156-06X	Pipette tip, Combo, 200ul , 6 stack (96 each), 576 tips

Note:

1. Do not touch pipette tips with bare hands to avoid RNase contamination.
2. Pipette tip box should be properly aligned on the instrument. Improper alignment may lead to display of “Error” messages.