


DATA SHEET
10ml Empty Assay Vial (Brown), 24/pk

Doc. No. 932- XT158-24X, Rev. No. B
Date of Release: 29-Apr-2021

Intended Use

This product is intended for laboratory use only. The 10ml Empty Assay Vial (Brown), 24/pk is for use on BioGenex *NanoVIP*[®] or *NanoMtrx*[®]100 Staining Systems.

Items Supplied:

Cat. No.	Description	Figure
XT158-24X	10ml Empty Assay Vial (Brown), 24/pk	

Note:

1. Proper positioning of all vials in reagent racks is crucial since the pipette head is calibrated to this height and orientation.
2. Recalibration of the reagent vial position is required when using the vials for the first time.
3. When loading these vials into the Automated Staining Systems reagent racks, make sure that the vial neck is to the left.