

**DATA SHEET****17ml Empty Assay Vial (Brown) for Xmatrix, 24/pk**

---


Doc. No. 932-XT160-24X Rev. A  
Date of Release: 19-Oct-2021

---

**Intended Use**

This product is intended for laboratory use only. The 17ml Empty Assay Vial (Brown), 24/pk is for use on BioGenex Xmatrix Infinity Staining Systems.

**Items Supplied:**

<b>Cat. No.</b>	<b>Description</b>	<b>Figure</b>
XT160-24X	17ml Empty Assay Vial (Brown), 24/pk	

**Note:**

1. Proper positioning of all vials in reagent racks is crucial since the pipette head is calibrated to this height and orientation.
2. Recalibration of the reagent vial position is required when using the vials for the first time.
3. When loading these vials into the Automated Staining Systems reagent racks, make sure that the vial neck is to the left.