

New

NanoMtrx[®]100

State-of-the-art
Fully Automated IHC and ISH System



NanoMtrx[®]100


NanoMtrx[®] 100 - compact 10 slide, fully automated benchtop user-programmable system for IHC and ISH

Introducing the revolutionary **micro-chambers[®]** module for uniform, reproducible, clean and intense stains



Key Benefits

- Multi-format specimen processing: FFPE, frozen, cell preparations, smears, and FNAs
- Micro-chamber® for uniform staining throughout the slide
- Temperature controlled micro-chambers® for minimal reagent consumption
- On-board auto-DAB mixing
- Gentle wash and blow-dry to eliminate tissue lift-off
- Generates 70% less waste
- Fast turnaround time of 2 hours with simultaneous 10 slide processing
- Intuitive user friendly GUI



NanoMtrx[®]100



Priming



Probe Wash



Slide Wash



Slide Blow




Micro-chamber Wash



Micro-chamber Blow Dry



Reagent Aspiration



Reagent Dispense

Contact Information:

USA: 48810 Kato Road, Suite 200E, Fremont, CA 94538, USA

India: APIIC SEZ, Adibatla, R.R.Dist. 501510, Telangana, India

China: Unit 511B, 5F Front Hall, Shanghai Exhibition Center,
Shanghai 200041, PRC, China

Customer service - US: customerservice@biogenex.com

Customer service - India: indiacs@biogenex.com

Customer service - Global: internationalcs@biogenex.com

Instrument support: instrument.care@biogenex.com

Reagent support: reagentsupport@biogenex.com

For US

T: +1-800-421-4149

F: +1-510-824-1490

For India

T: +91-9100032862

T: +91-6309448856

For China

T: +86-21-2220-6680

F: +86-21-2220-6679

For outside US, India, and China

T: +91-7702301043



www.biogenex.com